VARILITE® Technology

VARILITE'S® Air-Foam Floatation Technology combines the best characteristics of air and foam!

AIR: It is fluid and displaces under load providing pressure relief.

FOAM: It acts like a spring and compresses under load providing support and conformation.


Air-Foam Floatation – Key is immersion. The valve releases air to immerse the user in the foam. Each cell of the open-cell foam provides support and distributes the pressure. Increased immersion provides increased load distribution and greater stability. Optimal pressure distribution is achieved over the largest possible area. The combination of soft foam under the pelvis and firmer foam under the thighs creates continually adjusting forces that evens out the interface pressures.

What's needed from a cushion used on a wheelchair?

1. Pressure distribution
2. Postural support
3. Vibration reduction
4. Light weight
5. Low maintenance

VARILITE® cushions excel in all of these areas. The combination of soft foam (green) under the pelvis and firmer foam (blue) under the thighs creates continually adjusting forces (orange arrows) that evens out the interface pressures (red arrows). Air supports the load, while foam keeps the air where it's needed and preserves the cushion shape.

49 Bespoke Inc., distributes **VARILITE®** products in Canada and has a wealth of information on-line at www.varilite.ca. Please check out our products and give us a call for any further information or in-service.

t. 416-661-4499
tf. 866-827-4548
f. 416-661-8903
tf. 877-731-7745
info@varilite.ca

49BESPOKE.COM

Providing information and education through our relationship with Westech Health Care, Ltd.
info@westechhealth.com

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