

The ultimate mobility tool

Batec handbikes are a comprehensive solution to your mobility problems:

- ✓ Go from reduced to **augmented mobility**; from needing help to being **independent** and **proactive**.
- ✓ Get a **2 in 1**, without giving up your manual wheelchair.
- ✓ **Overcome obstacles** and **drive over bumpy terrains**.
- ✓ Attach our handbikes **in just seconds** and **with no need for transfers**.
- ✓ Also available for **quadriplegics** and **hemiplegics**.

You CAN get used to being in a **wheelchair**



But not to giving up on

LIFE

www.batec-mobility.com

49
BESPOKE

Distributed in Canada by
49 Bespoke Inc.
137 Buttermill Avenue,
Concord, Ontario L4K 3X5
www.49bespoke.com
Ph: +1.866.827.4548
Fax: +1.877.731.7745
Email: 49@49bespoke.com



PARA
QUAD
HEMI
version available

BATEC
MOBILITY



BATEC ELECTRIC 2

The perfect balance

BATEC ELECTRIC 2, the perfect balance

The **BATEC ELECTRIC** is the latest in our saga of handbikes. It features the latest **Batec Mobility** technology and an **18" Wheel** designed exclusively for this model, which means you can use it both in **urban** and **off-road settings**.



All-in-one LCD screen with gear selector.



18" Wheel, ultra-grip tyre and aluminium double-wall rim.



250 RPM 900W high torque brushless motor with reverse gear. Maximum speed 25 Km/h.



180 mm disc, two mechanical brakes, two aluminium levers with parking brake.



700-lumen front light with position lights / high beam options. Dual LED tall light on the stand frame.



Horn and USB port.



Removable right dropout for easy inner tube or tyre change.



Removable QR weights.



21,5 km* autonomy with 36 V – 280 Wh battery. 44 km* autonomy with 36 V – 576 Wh battery. 4A charger.



Matte black aluminium handlebar with folding system for convenience.




G&B handlebar option without levers for quadriplegics.



Handlebar option with right or left Monolever brakes for hemiplegics.

MEASUREMENTS



COLOURS



* Autonomy calculated under the following conditions: 75 - kg user, ultra-light wheelchair, flat terrain, no wind and steady speed of 15 km/h. May vary +/- 20% depending on the battery.
Number/version: MM-008/8.0. | July 2024